

Technical Specifications

Chambers	3
Compression Mode	4
Time Range	1min-24hours or continuous
Pressure Range	40-150mmHg
Hold Time	0-10s(Step by 1s)
Inteval Time	10-60s(Step by 1s)
Interface	5" touch screen
Rechargeable Battery	Yes (Lithium battery : 11.1V 9000mAh) Maximum Running Time: No Less than 8 Hours
Adaper	Yes
Input Voltage	AC100-240V 50-60Hz
Noise Level	<50 dB
Dimension	250(L)x184(W)x150(H)mm
Weight	2.5kg

LGT-2200DVT

Compression Therapy Device



Longest
PHYSIO & REHAB DEVICES

Guangzhou Longest Medical Technology Co., Ltd.

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999 **Fax:** 86-020-66353920

E-mail: export@longest.cn **Website:** www.longestmedical.com



Health Life Longest Care


Longest



LGT-2200DVT is a special Pressotherapy system that is designed to decrease the risk of DVT(Deep Venous Thrombosis) and Pulmonary Embolism(PE)before, during and after surgery


4 compression modes with reusable & disposable sleeves meet all your need

⊕ Compression Modes



⊕ Introduction of deep vein thrombosis (DVT)

⊙ Reason




Deep vein thrombosis is one of the most common complications of bed rest, prevention is urgently needed

⊕ Features

- ⊙ 4 Compression Modes
- ⊙ 5" Touch Screen
- ⊙ Reusable & Disposable Sleeves
- ⊙ Gradient Pressure Setting
- ⊙ Intermittent Pneumatic Compression
- ⊙ Pressure of Every Chamber Adjustable
- ⊙ Rechargeable Battery
- ⊙ Powered by AC or DC



⊙ For Disposable Sleeve



*Notice: For one - leg treatment, please use single - head hose.

⊙ For Reusable Sleeve



*Notice: For two - leg treatment, please use double - head hose.