

Certification



State Food and Drug Administration



CE (MDD)



ISO13485 / ISO9001



NRTL

Sleeves

Picture	Sleeve Type	Order NO.	Chamber
	Arm Sleeve	220041	4
		220081	8
	Leg Sleeve	220042	4
		220082	8
	Waist Sleeve	220083	4
		220083L	
	Full Body Sleeve	220016FB	8



LGT-2200WM

Pressotherapy Device for Medical Aesthetics



Guangzhou Longest Medical Technology Co., Ltd.

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999

Fax: 86-020-66353920

E-mail: export@longest.cn

Website: www.longestmedical.com



Health Life Longest Care


Longest



LGT-2200WM


Pressotherapy, For Your Future Beauty

LGT-2200WM is a sequential pneumatic pressure system. It realises therapeutic effects through program-controlled inflation and deflation of the sleeves covered in arms, legs, and waist.



Specifications

- ⊗ Channel: 4
- ⊗ Chambers: 4/8
- ⊗ Compression modes: 8
- ⊗ Time range: 1-99 min
- ⊗ Pressure range: 37-187 mmHg
- ⊗ Interface: 5.7" Touch screen
- ⊗ Weight: 10 kg
- ⊗ Dimension: 366mm(W) x 393mm(L) x 194mm(H)



Mechanism & Functions

The multi-lymphatic drainage system LGT-2200WM is used to drain the dissolved or melted fat cells generated after various obesity equipment by improving lymph glands or blood vessels that are constricted or damaged due to cellulite or excessive fat cells on the lymphatic drainage step, one of the general obesity lipolysis steps.



Features

- ⊗ User-friendly interface, just a few steps to finish your setting
- ⊗ Gentle inflation and deflation
- ⊗ Durable & easy-use sleeves with overlapping chambers
- ⊗ Pressure unit switchable between mmHg and kPa
- ⊗ 8 compression modes, 4 output channels
- ⊗ Real-time pressure, and inflation part monitoring
- ⊗ 2 pcs of 8-chamber cuffs can be connected with one device at the same time

Application

- ⊗ Lymph drainage
- ⊗ Blood circulation
- ⊗ Detoxification
- ⊗ Cellulite treatment
- ⊗ Body shaping in legs, abdomen, buttocks and arms
- ⊗ Strengthening of loose skin
- ⊗ Relaxation of stress
- ⊗ Relieve fatigue